


Face mask PP n-w 3 ply with elastic blue

Artikel number :	DK03B50AC
Quality :	EN 14683:2019+AC : 2019
Colour :	Blue
Inner material :	Nonwoven polypropylene
Kleur :	blue
Weight of Material :	22 gsm \pm 1%
Middle material :	Meltblown non-woven
Colour :	White
Weight of Material :	25 gsm \pm 1%
Other material :	Elastic band
Colour :	White
Material :	Elasthane



Dimensions of finished product

Length :	23 cm \pm 0,5cm
Height :	10 cm \pm 0,5cm
Packaging :	48 boxes with 50 pieces.

Method of use : Press the nose bar on both sides of nose to make the shape of the bar fit the nose. Fasten the loop on ears and pull the down edge of face mask to cover chin of user.

Warning: To avoid danger of suffocation, keep this wrapper away from babies and children.

10/2020